


M

Automatic Solvent Recognition wizard: example of an AI at work

Stan Sykora, Carlos Cobas,
Felipe Seoane, Pablo Monje, Esther Vaz, et al




M

GSD – Global Spectral Deconvolution

f-domain algorithm which automatically decomposes sets of superposed *Lorentzians* and *near-Lorentzians* and ends up in a

GSD Peaks List

GSD is a standard feature of Mnova,
subject to continuous development



GSD has been born in 2008

and presented in detail at several meetings in 2009 and 2010



MESTRELAB RESEARCH
Chemistry Software Solutions

M

GSD Peaks Editing

Peaks editing brings in
specific NMR know-how
by **classifying** of the peaks in an GSD List

GSD Peaks

Report Copy Set Flags Select Peaks Reveal Peaks Setup

Resolution: 1.00 [Normal]; Refinement: 2; AutoEdit: yes

	ppm	Width	Height	Area	Type	Flags
5	9.247	0.972	7.571	99.762	Compound	Weak
6	9.226	0.866	946.878	10952....	Compo...	None
7	9.222	1.004	1151.325	13926....	Compo...	None
8	9.217	0.950	1018.984	11608....	Compo...	None
9	9.214	1.287	1059.235	17003....	Compo...	None
10	9.196	4.133	25.819	1114.023	Impurity	Weak + Labile
11	8.850	1.461	6.746	111.421	Impurity	Weak





GSD Peaks Editing Modes

- ✓ Manual: Always accessible
- ✓ Automatic: An optional GSD feature always used in ASV

The screenshot displays the GSD Peaks software interface. The main window shows a table of peaks with columns for ppm, Width, Height, Area, and Flags. The 'GSD Settings' dialog box is open, showing various refinement options. The 'Auto edit' checkbox is checked and highlighted with a red box. The 'GSD Peak Flags' dialog box is also open, showing a dropdown menu for 'Type' with 'Solvent' selected.

GSD Settings

Settings

Refinement Level
Ref. 1 (2 fitting cycles)

Resolution

High
 Normal
 Low
 Custom 1.00

Additional Refinement

Splittings

Auto edit

Set as Default
Restore

GSD Peak Flags

Peak: 9.214

Type: Solvent

Flags: S. Reference, Q. Reference

OK
Cancel

ppm	Width	Height	Area	Flags
5	9.247			Peak
6	9.226			one
7	9.222			one
8	9.217			one
9	9.214	1.287	003....	Compo... None

AutoEdit does not exclude subsequent manual peaks editing
(transparently followed by **PostEdit**)



The GSD Auto-Edit algorithm: what is it up to

- ✓ Reference peaks recognition
- ✓ ^{13}C satellites recognition (where possible)
- ✓ Solvent peaks recognition (primary and secondary)
- ✓ Labile peaks recognition
- ✓ J-coupling multiplets recognition
- ✓ Impurity peaks recognition
- ✓ Weak peaks labeling





GSD Auto Editing results

Peak Types:

pt_Compound,
pt_Artifact,
pt_Impurity,
pt_Solvent,
pt_SReference,
pt_QReference

Peak Flags:


pf_SetByUser
pf_Hidden
pf_Weak
pf_C13Satellite
pf_Rotational
pf_Labile
pf_HasC13Sats





What can one do with the Edited Peaks List


a) Color-labeled peak plots of the edited peaks






What can one do with the Edited Peaks List

b) Synthetic spectra of any desired category of peaks





What can one do with the Edited Peaks List ?

c) Good editing is an essential pre-requisite of ...

... **ASV !**

Automatic Structure Verification



Solvent recognition

Possibly the most important and tricky part of **AutoEdit**

Human solvent recognition:

Based on a few imperfect rules and lots of intuition

Automatic solvent recognition:

Based on a long list of imperfect scoring tests

Primary and secondary solvent signals:

Primary solvent recognition is relatively easy

Secondary solvent signals (water) often present problems

Basic development rule:

How the hell does Tony know that this peak belongs to solvent !!!

(if you ask him, he does not know)





Solvent description

```
DWORD SdFlags; // See the FFP_Flags enum in EbPeaks.h
// Main Multiplet:
REAL Shift; // typical main multiplet shift value
REAL ShiftLoPpm; // low limit of main multiplet in ppm
REAL ShiftHiPpm; // high limit of main multiplet in ppm
DWORD Npeaks; // number of main peaks (may not be 0)
REAL JHz; // coupling constant
REAL Pweights[5]; // array of relative weights of the peaks
// Satellite Peaks:
REAL SatWeight; // satelites relative weight (0.011 or 0)
REAL JSatHz; // satelites coupling constant
REAL SatShiftPpm; // satelites relative isotope shift
// Associated Multiplet:
DWORD Npeaks2; // number of peaks (may bee 0)
REAL Weight2; // relative weight (typical)
REAL ShiftPpm2; // relative isotope shift
REAL JHz2; // coupling constant
REAL Pweights2[5]; // array of relative weights of the peaks
```





Solvent descriptors

Typically, a solvent entry in Mnova contains two Solvent Descriptors:

- One for the **primary solvent** peaks and
- One for the **secondary solvent** peaks (usually water)

However, it could contain any number of them (even a typical impurity might be nominally considered part of the solvent)

Coming soon:

Possibility to easily edit and add solvent description XML files.






Solvent recognition algorithm

- ✓ Look for the central peak of a solvent multiplet
- ✓ Scan all peaks in the admissible range
- ✓ Associate a full-fledged scoring system with each peak
- ✓ Scan on all of these criteria (if pertinent):
 - Position
 - Width
 - Height
 - Kurtosis
 - Presence of J-coupling multiplet peaks
 - J-HD coupling constant value
 - Relative intensities
 - Presence of ^{13}C satellite peaks
 - Relative intensities
 - J-CH coupling value
 - J-CH isotope shift
 - Presence of secondary J-coupling multiplet
 - Relative intensities
 - J-HD coupling in secondary multiplet




M

Some examples



M


Some examples



MESTRELAB RESEARCH
Chemistry Software Solutions



M



Thank You for your Patience

Any Questions ?



MESTRELAB RESEARCH
Chemistry Software Solutions

www.mestrelab.com

SMASH 2010, Portland, OR, USA

DOI: [10.3247/SL3nmr10.007](https://doi.org/10.3247/SL3nmr10.007)

www.ebyte.it

